ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE: The intended user should be a licensed medical professional, the patient, the patient's caretaker, or a family member providing assistance. The user should be able to understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

INTENDED USE / INDICATIONS: The Air-Stirrup Plus™ provides support to prevent further injury by effectively stabilizing the ankle. Soft-goods/semi-rigid designed to restrict motion through elastic or semi-rigid construction.

CONTRAINDICATIONS: N/A

WARNINGS AND PRECAUTIONS: 1) Do not use over open wounds. 2) Do not use if you experience any sensation or other unusual reactions occur while using this product. You should contact your doctor immediately. 3) Do not use this device if it was damaged and/or package has been opened. 4) If you develop an allergic reaction and/or experience itchy, red skin after coming into contact with any part of this device, please stop using it and contact your healthcare professional immediately.

NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device. If you have any of the following signs and symptoms, contact your healthcare professional immediately: • Numbness or tingling in your foot. • Cold or pale toes. (If you are wearing the Ankle Wrap, loosen or remove it. If the problem persists, contact your healthcare professional.) • Inability to put weight on your injured foot four days after injury. • Excessive pain that does not improve with rest. • Increase in swelling along with a decrease in mobility.

STEP 1: REDUCE PAIN AND SWELLING

1) Raise your foot. 2) Wrap cold pack in a moist cloth and place on your ankle. (Fig. 1) 3) Use 15 minutes at a time. The Ankle Wrap (included in the Air-Stirrup Plus) will help reduce swelling and should be used for the first few days only. The wrap does not provide support.

STEP 2: APPLY THE AIR-STIRRUP ANKLE BRACE, THEN WALK

1) Open brace and lift bottom of arches until heel straps are fully exposed. 2) Peel up heel straps and adjust heel pad width for a snug fit. (Fig. 3) 3) Firmly press heel straps and arches back in place. 4) Secure straps from bottom to top. 5) Place round edge of heel pad under heel. 6) Align brace sides with ankle. 7) Secure straps from bottom to top. 8) Put on lace shoe or sneaker.

ANKLE BRACE FIT ADJUSTMENT

If you experience pinching or uneven pressure, reposition brace or readjust heel pad and straps. Ankle brace arches are preinflated and normally do not require adjustment.

HIGH ALTITUDE: When flying, readjust the straps to a comfortable pressure. At high altitudes the ankle will expand beyond their optimal level and air will have to be removed from the air cells. Insert inflation tube about 1/2" into valve on top of arches. Remove a small amount of air by gently squeezing arches, or a small amount of air blowing into tube (from mouth) as you release pressure. Seal valve by squirting the valve flat. Wait below end, and slowly remove tube. Fold and close valve.

Correctly position your ankle. Your injured ankle must be correctly positioned with the brace on, before standing or walking. While seated, bend your knee and slowly move your injured foot in toward your body until your foot is flat on the floor.

Step 3: Exercise

When you begin these exercises, you may have additional soreness — this is normal. Progress through the exercises as your ankle will allow. If pain becomes excessive — do not subside with rest and ice, contact your healthcare professional immediately (see "Warning Signs" on inside back cover).

Toe Point Extend your leg and slowly point your toes forward, then slowly point your toes toward ceiling. (Fig. 3)

Alphabet Exercise Write your name on your injured leg on top. Draw the alphabet in capital letters with your big toe. (Fig. 3)

Balance Exercises Motion Exercises should be started immediately.

Achilles Stretch Each exercise should be done 10 times and repeated 3 times a day, or as recommended by your healthcare professional.

Balance Exercises After you have successfully completed the Motion Exercises, begin the Balance Exercises. For each exercise you will have to stand, using a chair or wall for balance. Do each exercise 10 times and repeat 3 times a day.

Toe Raise Stand on both feet, holding a chair or wall for balance. Slowly move up and down on the toes of both feet. Raise your uninjured foot. Slowly move up and down on the toes of your injured foot. (Fig. 4)

Single Leg Stand Stand on your injured foot flat on the floor and your other foot raised. Balance for up to 30 seconds.

Care The entire Air-Stirrup Ankle Brace can be hand washed in lukewarm water with mild soap and air dried.

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

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NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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